

BAILEY

Model 4311 Cabinet Treatment Table with Adjustable Dual Legrest & Adjustable Backrest

Model 4311 Treatment Table Features:

- 30" height, 30" wide and 78" long
- 2" thick sewn vinyl upholstery in US 387 Camel standard, with your choice of Bailey upholstery colors, and either Natural or Dark Walnut wood finishes (no cost), or optional Sherwin Williams colors (upcharge).
- Four Doors with Wire Fingerpulls (Classic Style), Optional Pro Style Doors with Recessed Pulls
- Heavy Duty Adjustable Legrests and Backrest, each with 8 Locking Positions (up to 65 degrees for Backrest, and 32 Degrees for Legrests)
- Weight capacity 350 lb.
- Ships fully assembled; ship weight 260 lb.



Model 4311 Cabinet Table
Features 4 Door Cabinet,
484 Legrests and 485 Backrest
2" Sewn Top shown in American Beauty
Upholstery, in Natural Wood Finish)

CORPORATE MEMBER

