

BAILEY

Custom Dual Whirlpool Tables

Dual Whirlpool Table Features:

- Choice of ANY tank width, up to 26" wide.
- Specify by actual tank **floor-to-rim** height in inches, 22" to 36" high tank.
- Upholstered Vinyl Top and Twin Combo Seats.
- Choice of Natural, Dark Walnut, or custom wood finishes.
- May mix two different WIDTH whirlpools, but must be SAME HEIGHT whirlpools.
- Top size: 32" Deep x 72" Wide minimum, 84" maximum (depends on tank sizes).
- Adjustable Height Foot Glides.
- Weight capacity 350 lb.
- Ship weight 250 lb. [330 lb. with shelves].
- Optional Dual Rear Shelves add storage.
- Upholstered 1" Padded Top protects against water damage.
- Other Custom variations available, ask for a quote!

Dual Whirlpool Table

[Shown in Natural, with Padded Top Option]



OPTIONS: Rear Shelves and Upholstered Top

[Shown in Royal Blue Vinyl and Natural Finish]



Corporate Member



**See our full catalog of products at
www.baileymfg.com**