

# **SCHOOL FOCUS PRODUCTS**

#### CLASSROOM CHAIR

The seat and foot board of this versatile chair can be adjusted to accommodate children of various sizes. This unit, styled to resemble a typical classroom chair, has a padded upholstered seat and sling style back, with skis to prevent tipping.

#### **MODEL 154 - ADJUSTABLE CHAIR**

- seat 12" wide, 14" from floor
- seat depth adjusts from 6" to 12"
- arm height 18-1/2" from floor
- footboard adjusts 5" to 11" from seat
- 58" body or leg strap with velcro
- chrome plated frame
- chestnut vinyl upholstery

#### MODEL 123 - ACCESSORY TRAY (SHOWN IN IMAGE OF MODEL 154)

- fits MODEL 154
- chrome plated frame
- chestnut vinyl upholstery





MODEL 154



**MODEL 1719** 

## MULTI-USE CLASSROOM CHAIR

Adjustable. Adaptable. Affordable. With the addition of easily attached accessories, this basic chair can be customized to meet the individual needs of any neurologically delayed or orthopedically involved child up to 48" tall. Standard color is chestnut.

#### **MODEL 1700 - CHILD SIZE**

- back can be raised and lowered
- footrest adjusts 5" to 11" from seat
- seat depth adjusts from 6" to 12"
- seat 11-1/2" wide x 14" deep
- seat belt included
- stabilizing skis prevent tipping
- 14" to seat. 18" to armrest

#### **MODEL 1701 - ADOLESCENT SIZE**

- back can be raised and lowered
- footrest adjusts 3" to 13" from seat
- seat depth adjusts from 9" to 15"
- seat 14" wide x 16" deep
- seat belt included
- stabilizing skis prevent tipping
- 17" to seat, 24" to armrest



MODEL 1723 & 1724



**MODEL 1727** 

### **ACCESSORIES**

MODEL 1719 - ADJUSTABLE LATERAL SUPPORTS- FOR MODEL 1700 AND MODEL 1701

MODEL 1727 - ADJUSTABLE KNEE ABDUCTOR- FOR MODEL 1700 AND MODEL 1701

MODEL 1723 - ADJUSTABLE PELVIC SUPPORTS- FOR MODEL 1700

**MODEL 1724 - ADJUSTABLE PELVIC SUPPORTS- FOR MODEL 1701** 

MODEL 1715 - CLEAR ACRYLIC TRAY- FOR MODEL 1700 (23" WIDE X 18" DEEP)

MODEL 1716 - CLEAR ACRYLIC TRAY- FOR MODEL 1701 (23" WIDE X 18" DEEP)

**MODEL 1731 - ADJUSTABLE FOOT BRACKET- FOR MODEL 1700** 

MODEL 1732 - ADJUSTABLE FOOT BRACKET- FOR MODEL 1701



MODEL 1731 & 1732



MODEL 1715 & 1716

#### BAILEY MANUFACTURING CO.

WWW.BAILEYMFG.COM



**Bailey Manufacturing Company** 

bailey\_mfg\_lodi

P.O. BOX 130 LODI, OHIO 44254-1030 IN USA: PH 800-321-8372 FAX: 800-224-5390 OUTSIDE USA: PH 330-948-1080 FAX 330-948-4439

Over 60 years of experience in crafting equipment for physical therapy, occupational therapy, work hardening, and sports medicine

#### **CONTACT:**