

# BAILEY

## Model 2420 Tee Stool

### **Features:**

- Five Adjustable Height Range positions from 12" to 16"
- Heavy Duty Rubber Foot will not mark floors
- 200 LB weight capacity
- Birch Poplar Hardwood Seat, 10 3/4" Diameter
- Durable Chrome Plated Steel Construction
- Encourages balance and develops integrative and perceptual motor skills.
- Suitable for use in special needs classrooms and for home use.
- Ships UPS; Assembly Required



**NEW**

**Heavy Duty Floor Tip  
(New Design Beginning Jan 2017)**

**Retrofit Kit (p/n A2420HD)  
is Available to upgrade our  
older Tee Stools**