

BAILEY


Manufacturing Company

The Quality Leader in Sports Medicine Products





Danny Poole - Director of Sports Medicine - Head Athletic Trainer at Clemson University



"The custom equipment Bailey helped us design into our West End Zone Athletic Training Room is top notch, and makes our facility very efficient in dealing with our Athletes."



Treatment Table Dimensions and Models

Table Styles	Upholstered Top 24" Wide 78" Long 1" Top	Upholstered Top 30" Wide 78" long 1" Top	Upholstered Top 30" Wide 78" long 2" Top	Raised Rim Top 30" Wide 78" long Acc. Mat optional
Plain Table	400	410	430	460
Table with H-Brace	401	411	431	461
Table with plain shelf	402	412	432	462
Table with divided shelf	403	413	433	463
Table with drawer and storage shelf	404	414	434	464
Table with drawer and H-Brace	405	415	435	465
Table with drawer and plane shelf	406	416	436	466
Table with drawer and divided shelf	407	417	437	467
Table with adjustable back		418	438	468
Table with adjustable back and drawer		419	439	469
Table with adjustable back and H-Brace		420	440	470
Table with adjustable back and plain shelf		421	441	471
Table with adjustable back and divided shelf		422	442	472
Table with adjustable back, H-Brace and drawer		423	443	473
Table with adjustable back, plain shelf and drawer		424	444	474
Table with adjustable back, divided shelf and drawer		425	445	475
Table with adjustable back, enclosed cabinet and sliding doors		426	446	476
Table with enclosed cabinet, sliding door and drawer		427	447	477



Model 432



Model 441



Model 434 with optional Model 40 plain shelf

Model 720

Model 750

Economy Hi-Low 2" Upholstered Top Model 4243

Upholstered shelf
30" x 78"
Black enameled steel frame
Adjustable back
Height adjustment from 25" - 33"
Manual hydraulic lifting system

Professional Hi-Low Treatment Table Model 4050

Adjusts from 27" to 39" high
Foot control
UL hospital grade plug
Natural finish is standard, Walnut finish is available.
30" x 78" x 2"
110V 60Hz 6.5AMP 720WATT Thermal overload protection



Model 4090

Professional Hi-Low Upholstered Mat Tables - Electric

All the durability and quality features that a professional demands for an efficient and economical department. A Foot Switch is standard. Electric height adjustment from 23" to 32" high. 110V 60Hz, UL hospital grade plug. 500 LB weight capacity
Model 9620 - 4' x 7'
Model 9621 - 5' x 7'
Model 9622 - 6' x 8'



Model 9620

Economy Electric Hi-Low 2" Upholstered Top Model 4250

30" x 78" top with
Full size Drawer and Gas Spring Adjustable Back with dual controls



Model 4050

Professional Hi-Low Treatment Table Model 4090

Adjusts from 27" to 39" high
Foot control
UL hospital grade plug
Natural finish is standard, Walnut finish is available.
Functional shelves
Adjustable backrest
30" x 78" x 2"
110V 60Hz 6.5AMP 720WATT Thermal overload protection

Upholstered Top Mat Tables - Electric Hydraulic

Height adjusts from 19" to 27" high
Black enamel steel frame
Built in floor levelers
Natural finish is standard
Walnut finish is available
Toggle Switch is Standard, a foot switch is optional
Model 4124 - 3' x 7' x 2"
Model 4125 - 4' x 7' x 2"
Model 4126 - 5' x 7' x 2"
Model 4127 - 6' x 8' x 2"



Model 4125

BariMatic Electric Hi-Lo Upholstered Mat Tables

20" to 28" height adjust, using Foot Switch (standard).
Adjustable floor levelers.
Camel color Soft Vinyl top cover standard, with optional Herculite and soft vinyl
All steel, heavy duty powder coated frame with crossbrace for stability.
Lifting Weight capacity 1000 lb.
Model 4554 - 4' x 7'
Model 4555 - 5' x 7'
Model 4556 - 6' x 8'



Model 4554



Model 4305

**Cabinet Table
Model 4305**

30" wide x 78" long x 1" top.
Six drawers.
One door.
Adjustable back.
Natural finish is standard.
Walnut finish is available.

**Cabinet Table
Model 4301**

30" wide x 78" long x 1" top.
Four doors.
Adjustable shelf.
Natural finish is standard.
Walnut finish is available.



Model 4301



**Model 7720 Dual
Shown with optional shelf**

**Whirlpool Tables
Models 7720 / 7724**

Choice of 20" wide tank [Model 7720] or 24"/25" wide tank [Model 7724].
Specify by actual tank floor-to-rim height in inches, 22" to 36" tank.
Upholstered Vinyl Top and Twin Combo Seats.
Choice of Natural or Dark Walnut wood finishes.
Available for either 20" or 24"/25" wide whirlpools.
Top size: 32"D x 48"W; Reinforced Combo Seats: 14" square.
Adjustable Height Foot Glides.
Weight capacity 350 lb.
Ship weight 125 lb. [170 lb. with optional shelves].
Optional Dual Rear Shelves add storage.
Upholstered 1" Padded Top protects against water damage.
Other Custom variations available, ask for a quote!
Available in Single, Dual (pictured) and Triple tables
Please specify actual floor to rim height in inches.
22"- 36" tanks accepted



Model 7724



**Need Whirlpool Chairs?
see page 15 for more details**

Adjustable Back Bariatric Treatment Table with H- Brace Model 4515 / 4516

2" upholstered top.
 Adjustable back features positive locking in 6-8 places
 Standard H-Brace
 Natural finish is standard.
 Walnut finish is available.
 6 Legs and special frame design allows 1000 pound patient capacity
 Version with Plain Shelf is available as Model 4516
 34" wide x 78" long x 28" high



Model 4515

Bariatric Mat Table Model 4520

6 Leg and special frame design allows 1000 pound patient capacity
 4' wide x 7' long x 18" high with 2" upholstered top
 US 387 Camel, standard.
 Choice of optional Bailey upholstery Colors
 Walnut finish is available.



Model 4520

Bariatric Treatment Table with Plain Shelf Model 4510 /4511

6 Leg and special frame design allows 1000 pound patient capacity
 34" wide x 78" long x 28" high
 2" thick upholstery top, camel standard
 Version with H-Brace available only as Model 4510
 Choice of optional Bailey upholstery Colors



Model 4511

STANDARD FABRIC COLOR SELECTOR - SOFT VINYL





Model 484

Leg rests adjust to 32 degrees, with 8 positions in 3" increments.

**Double Adjustable Leg Rest Table
Model 484**

2" upholstered top.
Natural finish is standard.
Walnut finish is available.
30" wide x 78" long x 30" high

**Professional Adjustable Back Treatment
Table
Model 4505**

Adjustable from either side, the back locks in any position from horizontal to 60°
31" total height with non-skid rubber feet
30" x 78" top with 2" upholstered top.
All steel, powder coated frame with storage shelf standard
Available Options: Model 39 Nose Slot, Model 46 Paper Holder and Model 48 Adjustable Paper Cutter
Weight capacity 500lb./227Kg



Model 4505

**Professional Treatment Table
Model 4500**

31" total height with non-skid rubber feet
30" x 78" top with 2" upholstered top.
All steel, powder coated frame with storage shelf standard
Available Options: Model 39 Nose Slot, Model 46 Paper Holder and Model 48 Adjustable Paper Cutter
Weight capacity 500lb./227Kg



Model 4500



Model 485

Adjustable Back Rest and bKnee Gatch Table Model 485

2" upholstered top.
Natural finish is standard.
Walnut finish is available.
30" wide x 78" long x 30" high

Heavy Duty Adjustable Backrest inclines 60 degrees, with 8 positions in 3" increments

Heavy Duty Table Model 350

Commonly used for hip extension exercises, hamstring curls, quadricep strengthening exercises, as well as upper extremity exercises.

Thigh support strap (removable)

Two permanently mounted chrome plated torque units with sealed ball bearings

Each torque unit has one adjustable 6" weight-holder and one 8" screw-in weight-holder

30" wide x 32" deep x 30" high

Top upholstered in vinyl

37.5 pounds per torque unit capacity



Model 350

Professional Exercise Table Model 345

Commonly used for hip extension exercises, hamstring curls, quadricep strengthening exercises, as well as upper extremity exercises. The tapered backrest is operated by gas spring and adjusts from 0° to 90°.

Adjustable headrest (removable)

Contoured handholds with non-slip grips

Thigh support strap (removable)

Two permanently mounted chrome plated torque units with sealed ball bearings

Each torque unit has one adjustable 6" weight-holder and one 8" screw-in weight-holder

Back height 27-1/2" to top of headrest

Top upholstered in Black vinyl

30" wide x 22" deep x 30" high

Back height 27-1/2" to top of headrest

Top upholstered in Black vinyl

37.5 pounds per torque unit capacity



Model 345

Models 4310 & 4311 Full Cabinet Style Sports Medicine Treatment Tables

30" height, 30" wide and 78" long

2" thick sewn vinyl upholstery in US 387 Camel standard, with your choice of Bailey upholstery colors, and either Natural or Dark Walnut wood finishes (no cost), or optional Sherwin Williams colors (upcharge).

Twin Sliding Doors on hardwood tracks with enclosed storage area (Model 4310)

Four Doors with Wire Fingerpulls (Classic) Style or Optional Recessed Pulls on ProStyle Doors (Model 4311)

Heavy Duty Adjustable Legrests and Backrest, each with 8 Locking Positions (up to 65 degrees for Backrest, and 32 Degrees for Legrests) Weight capacity 350 lb.

Ships fully assembled; ship weight 260 lb.

Legrests may be ordered as LH or RH



Model 4311



Model 4310



Model 10

**Fully Enclosed Taping Table
Model 10**

Standard 1" Padded Upholstered Top.
Optional 1 1/4" wood grain laminate top. (shown)
3/4" multi-layered birch sides and shelves.
Standard height is 36"
Optional height of 30" or 42" available at no extra charge.
26" wide x 40" long.
Natural finish is standard. (Shown)
Walnut finish is available.
Upholstered top available as a special request item.

End Shelved Taping Table Model 12

Standard 1" Padded Upholstered Top.
Two End Shelves for storage.
Standard height is 36".
Optional height of 30" or 42" available at no extra charge.
Seat Top is 26" wide x 40" long.
Natural finish is standard.
Walnut finish is available. (Shown)
Custom sewn vinyl Team Color Stripe available as a low cost option. (Shown)
Custom Laminated Top available as an added cost option.



Model 12



Model 14


**Taping Table
Model 14**

Upholstered top 24" x 48" x 1".
Standard height is 36".
Optional height of 30" or 42" available at no extra charge.
Natural finish is standard. (Shown)
Walnut finish is available.
Consult the fabric color selector for the complete selection of colors available.

**Cuff Weights
and Weight
Wagons**
see page 18-
19 for more
details

**Deluxe Taping Table
Model 25**

Height choices: 30", 36" or 42"
1" thick top with vinyl upholstery, with your choice of optional Bailey upholstery colors
Natural or Dark Walnut wood finishes
Full Rear Shelf and End Storage Shelves
Stationary 60" Back support
Weight capacity 400lb
Shown with Custom Embroidery option




Model 25

Modular Taping Cabinets (note: use with Supply Cabinets to make “stations”)

Common Features: top 6” front overhang, 1” side overhang, two doors, adjustable inside shelf, backrest adds 12” to overall height

Model	Description	Height	Width	Depth
8001-C	Classic Taping Cabinet	36”	24”	36”
8002-C	Classic Taping Cabinet	42”	24”	36”
8003-C	Classic Taping Cabinet	36”	24”	42”
8004-C	Classic Taping Cabinet	42”	24”	42”
8005-C	Classic Taping Cabinet w/Backrest	36”	24”	42”
8006-C	Classic Taping Cabinet w/Backrest	42”	24”	42”

Model	Description	Height	Width	Depth
8001-P	Pro Taping Cabinet	36”	24”	36”
8002-P	Pro Taping Cabinet	42”	24”	36”
8003-P	Pro Taping Cabinet	36”	24”	42”
8004-P	Pro Taping Cabinet	42”	24”	42”
8005-P	Pro Taping Cabinet w/Backrest	36”	24”	42”
8006-P	Pro Taping Cabinet w/ Backrest	42”	24”	42”



Modular Supply Cabinets (note: use with Taping Cabinets to make “stations”)

Common Features: top 2” front overhang, sides flush, 3 drawers, “English Oak” laminate wastebin top

Model	Description	Height	Width	Depth
8021-C	Classic Supply Cabinet, 3 Drawers	30”	18”	36”
8022-C	Classic Supply Cabinet, 3 Drawers	30”	18”	42”
8023-C	Classic Supply Cabinet, 3 Drawers,	36”	18”	36”
8024-C	Classic Supply Cabinet, 3 Drawers,	36”	18”	42”
8025-C	Classic Supply Cabinet, 3 Drawers	30”	24”	36”
8026-C	Classic Supply Cabinet, 3 Drawers	30”	24”	42”
8027-C	Classic Supply Cabinet, 3 Drawers	36”	24”	36”
8028-C	Classic Supply Cabinet, 3 Drawers	36”	24”	42”

Model	Description	Height	Width	Depth
8021-P	Pro Supply Cabinet, 3 Drawers	30”	18”	36”
8022-P	Pro Supply Cabinet, 3 Drawers	30”	18”	42”
8023-P	Pro Supply Cabinet, 3 Drawers	36”	18”	36”
8024-P	Pro Supply Cabinet, 3 Drawers	36”	18”	42”
8025-P	Pro Supply Cabinet, 3 Drawers	30”	24”	36”
8026-P	Pro Supply Cabinet, 3 Drawers,	30”	24”	42”
8027-P	Pro Supply Cabinet, 3 Drawers	36”	24”	36”
8028-P	Pro Supply Cabinet, 3 Drawers	36”	24”	42”



Model	Description	Height	Width	Depth
8029-C	Classic Supply Cabinet, 1 Drawer, 1 Door	30”	18”	36”
8030-C	Classic Supply Cabinet, 1 Drawer, 1 Door	30”	18”	42”
8031-C	Classic Supply Cabinet, 1 Drawer, 1 Door	36”	18”	36”
8032-C	Classic Supply Cabinet, 1 Drawer, 1 Door	36”	18”	42”
8033-C	Classic Supply Cabinet, 1 Drawer, 2 Doors	30”	24”	36”
8034-C	Classic Supply Cabinet, 1 Drawer, 2 Doors	30”	24”	42”
8035-C	Classic Supply Cabinet, 1 Drawer, 2 Doors	36”	24”	36”
8036-C	Classic Supply Cabinet, 1 Drawer, 2 Doors	36”	24”	42”

Model	Description	Height	Width	Depth
8029-P	Pro Supply Cabinet, 1 Drawer, 1 Door	30”	18”	36”
8030-P	Pro Supply Cabinet, 1 Drawer, 1 Door	30”	18”	42”
8031-P	Pro Supply Cabinet, 1 Drawer, 1 Door	36”	18”	36”
8032-P	Pro Supply Cabinet, 1 Drawer, 1 Door	36”	18”	42”
8033-P	Pro Supply Cabinet, 1 Drawer, 2 Doors	30”	24”	36”
8034-P	Pro Supply Cabinet, 1 Drawer, 2 Doors	30”	24”	42”
8035-P	Pro Supply Cabinet, 1 Drawer, 2 Doors	36”	24”	36”
8036-P	Pro Supply Cabinet, 1 Drawer, 2 Doors	36”	24”	42”

Modular Taping Stations, Two Seats

(note: made with 8000 Series cabinets)

Common Features: top 6" front overhang, max; 1" side overhang

Model	Description	Height	Width	Depth
8201-SC	Classic Taping Station - 2 Seats	36"	66"	36"
8202-SC	Classic Taping Station - 2 Seats	42"	66"	36"
8203-SC	Classic Taping Station - 2 Seats	36"	66"	42"
8204-SC	Classic Taping Station - 2 Seats	42"	66"	42"
8205-SC	Classic Taping Station - 2 Seats	36"	72"	36"
8206-SC	Classic Taping Station - 2 Seats	42"	72"	36"
8207-SC	Classic Taping Station - 2 Seats	36"	72"	42"
8208-SC	Classic Taping Station - 2 Seats	42"	72"	42"


Model	Description	Height	Width	Depth
8209-SC	Classic Taping Station, w/Backrests - 2 Seats	36"	66"	42"
8210-SC	Classic Taping Station, w/Backrests - 2 Seats	42"	66"	42"
8211-SC	Classic Taping Station, w/Backrests - 2 Seats	36"	72"	42"
8212-SC	Classic Taping Station, w/Backrests - 2 Seats	42"	72"	42"

Model	Description	Height	Width	Depth
8209-SP	Pro Taping Station, w/Backrests - 2 Seats	36"	66"	42"
8210-SP	Pro Taping Station, w/Backrests - 2 Seats	42"	66"	42"
8211-SP	Pro Taping Station, w/Backrests - 2 Seats	36"	72"	42"
8212-SP	Pro Taping Station, w/Backrests - 2 Seats	42"	72"	42"


Model	Description	Height	Width	Depth
8201-SP	Pro Taping Station - 2 Seats	36"	66"	36"
8202-SP	Pro Taping Station - 2 Seats	42"	66"	36"
8203-SP	Pro Taping Station - 2 Seats	36"	66"	42"
8204-SP	Pro Taping Station - 2 Seats	42"	66"	42"
8205-SP	Pro Taping Station - 2 Seats	36"	72"	36"
8206-SP	Pro Taping Station - 2 Seats	42"	72"	36"
8207-SP	Pro Taping Station - 2 Seats	36"	72"	42"
8208-SP	Pro Taping Station - 2 Seats	42"	72"	42"
8209-SP	Pro Taping Station - 2 Seats	36"	66"	42"
8210-SP	Pro Taping Station - 2 Seats	42"	66"	42"
8211-SP	Pro Taping Station - 2 Seats	36"	72"	42"
8212-SP	Pro Taping Station - 2 Seats	42"	72"	42"



8210-SP with logo & custom stain



8312-SC with logo & Walnut stain



Classic Modular Taping Stations, Three Seats

(note: made with 8000 Series cabinets)


Common Features: top 6" front overhang, max; 1" side overhang

Model	Description	Height	Width	Depth
8301-SC	Classic Taping Station - 3 Seats	36"	108"	36"
8302-SC	Classic Taping Station - 3 Seats	42"	108"	36"
8303-SC	Classic Taping Station - 3 Seats	36"	108"	42"
8304-SC	Classic Taping Station - 3 Seats	42"	108"	42"
8305-SC	Classic Taping Station - 3 Seats	36"	120"	36"
8306-SC	Classic Taping Station - 3 Seats	42"	120"	36"
8307-SC	Classic Taping Station - 3 Seats	36"	120"	42"
8308-SC	Classic Taping Station - 3 Seats	42"	120"	42"

Model	Description	Height	Width	Depth
8309-SC	Classic Taping Station, w/Backrests - 3 Seats	36"	108"	42"
8310-SC	Classic Taping Station, w/Backrests - 3 Seats	42"	108"	42"
8311-SC	Classic Taping Station, w/Backrests - 3 Seats	36"	120"	42"
8312-SC	Classic Taping Station, w/Backrests - 3 Seats	42"	120"	42"

Model	Description	Height	Width	Depth
8301-SP	Pro Taping Station - 3 Seats	36"	108"	36"
8302-SP	Pro Taping Station - 3 Seats	42"	108"	36"
8303-SP	Pro Taping Station - 3 Seats	36"	108"	42"
8304-SP	Pro Taping Station - 3 Seats	42"	108"	42"
8305-SP	Pro Taping Station - 3 Seats	36"	120"	36"
8306-SP	Pro Taping Station - 3 Seats	42"	120"	36"
8307-SP	Pro Taping Station - 3 Seats	36"	120"	42"
8308-SP	Pro Taping Station - 3 Seats	42"	120"	42"

Model	Description	Height	Width	Depth
8309-SP	Pro Taping Station, w/Backrests - 3 Seats	36"	108"	42"
8310-SP	Pro Taping Station, w/Backrests - 3 Seats	42"	108"	42"
8311-SP	Pro Taping Station, w/Backrests - 3 Seats	36"	120"	42"
8312-SP	Pro Taping Station, w/Backrests - 3 Seats	42"	120"	42"



Classic Modular Taping Stations, Four Seats

Common Features: top 6" front overhang, max; 1" side overhang

Model	Description	Height	Width	Depth
8401-SC	Classic Taping Station - 4 Seats	36"	150"	36"
8402-SC	Classic Taping Station - 4 Seats	42"	150"	36"
8403-SC	Classic Taping Station - 4 Seats	36"	150"	42"
8404-SC	Classic Taping Station - 4 Seats	42"	150"	42"
8405-SC	Classic Taping Station - 4 Seats	36"	168"	36"
8406-SC	Classic Taping Station - 4 Seats	42"	168"	36"
8407-SC	Classic Taping Station - 4 Seats	36"	168"	42"
8408-SC	Classic Taping Station - 4 Seats	42"	168"	42"



8412 Pro with logo


Model	Description	Height	Width	Depth
8409-SC	Classic Taping Station, w/Backrests - 4 Seats	36"	150"	42"
8410-SC	Classic Taping Station, w/Backrests - 4 Seats	42"	150"	42"
8411-SC	Classic Taping Station, w/Backrests - 4 Seats	36"	168"	42"
8412-SC	Classic Taping Station, w/Backrests - 4 Seats	42"	168"	42"



Model	Description	Height	Width	Depth
8401-SP	Pro Taping Station - 4 Seats	36"	150"	36"
8402-SP	Pro Taping Station - 4 Seats	42"	150"	36"
8403-SP	Pro Taping Station - 4 Seats	36"	150"	42"
8404-SP	Pro Taping Station - 4bSeats	42"	150"	42"
8405-SP	Pro Taping Station - 4 Seats	36"	168"	36"
8406-SP	Pro Taping Station - 4 Seats	42"	168"	36"
8407-SP	Pro Taping Station - 4 Seats	36"	168"	42"
8408-SP	Pro Taping Station - 4 Seats	42"	168"	42"



Model	Description	Height	Width	Depth
8409-SP	Pro Taping Station, w/Backrests - 4 Seats	36"	150"	42"
8410-SP	Pro Taping Station, w/Backrests - 4 Seats	42"	150"	42"
8411-SP	Pro Taping Station, w/Backrests - 4 Seats	36"	168"	42"
8412-SP	Pro Taping Station, w/Backrests - 4 Seats	42"	168"	42"



Modality Cabinets & Carts

Model	Description	Height	Width	Depth
8601-C	Classic Modality Cabinet	30"	18"	20"
8602-C	Classic Modality Cabinet	32"	24"	20"
8603-C	Classic Modality Cabinet	30"	24"	20"
8604-C	Classic Modality Cart	30"	18"	20"
8605-C	Classic Modality Cart	32"	24"	20"
8606-C	Classic Modality Cart	30"	24"	20"

Model	Description	Height	Width	Depth
8601-P	Pro Modality Cabinet	30"	18"	20"
8602-P	Pro Modality Cabinet	32"	24"	20"
8603-P	Pro Modality Cabinet	30"	24"	20"
8604-P	Pro Modality Cart	30"	18"	20"
8605-P	Pro Modality Cart	32"	24"	20"
8606-P	Pro Modality Cart	30"	24"	20"



8603 Pro Modality Cabinet



NOTES:

- 1) All Modular cabinets and products are finished all around, and may be joined together or used freestanding; all cabinets have 4 leveling glides.
- 2) Taping Cabinet Table with Backrest adds 6" to depth and 12" to total height
- 3) Supply Top Wastebins are 28 quart, color is black
- 4) Natural Finish is standard; Walnut or other stained Finishes available at added cost. Stain color specified to Sherwin-Williams Number, or custom matched by request.
- 5) Tops: All standard Bailey Upholstery Colors are available, or can be special ordered. "English Oak" is the Standard laminate, other laminate colors are available, and must be specified by Brand and or level of quality.
- 6) Dimensions above are floor FOOTPRINT. There will be top overhangs of 1" each side of taping, 2" front of supply, and 2" to 3" front of taping cabinets.

Custom Options (most popular items)

Model	Description
8900	Lock Option - per Door
8901	Lock Option - per Drawer
8902-484	Adds Locking HD Dual Legrests to Treatment Table
8903-485	Adds Locking HD Adjustable Backrest to Treatment Table
8904-49	Adds Sliding Door Cabinet Base to Treatment Table
8905	Pressed 2" Thick Pad Upgrade - For 1" Taping or Treatment Top
8906	Sub a Top Tray for Top Drawer in Supply Cabinet
8907	Classic Fronts with Fingerpull Edges (No Hardware)
8908	Sherwin Williams Wood Finish (hand applied stain)
8909	Embroidered Logo in Backrest, Pillow or Table Top
8910	Custom Supply Top Laminate Color (specify by number)
8911	Taping Top Laminated (front only with backrest)
8912	Classic Front - Sub Brass Pulls for Chrome
8913	Sewn Multicolor Stripe in Upholstery (two colors)
8914	Standard Dark Walnut Wood Finish (machine applied)
8915	1" THK PAD, 26W x 30H, WALL BKTS INCLUDED



Model 26 With Option 8913 Multicolor Strip in Upholstery

Modular Taping Cabinets Side Storage Options

Model	Description
8916	Left Side Classic Taping Cab Side Storage - Full 4 Drawers, 1 Door
8917	Right Side Classic Taping Cab Side Storage - Full 4 Drawers, 1 Door
8918	Left Side Pro Taping Cab Side Storage - Full 4 Drawers, 1 Door
8919	Right Side Pro Taping Cab Side Storage - Full 4 Drawers, 1 Door
8920	Left Side Taping Cabinet Side Storage - Back Half, Open Shelf
8921	Right Side Taping Cabinet Side Storage - Back Half, Open Shelf
8922	Left Side Classic Taping Cabinet Side Storage - Back Half, 1 Door
8923	Right Side Classic Taping Cabinet Side Storage - Back Half, 1 Door
8924	Left Side Pro Taping Cabinet Side Storage - Back Half, 1 Door
8925	Right Side Pro Taping Cabinet Side Storage - Back Half, 1 Door
8926	Supply Cabinet, sub Full Extension Pull Out Tray for Top Drawer

Specify Left or Right Side as Facing Cabinet



8916 Pro

Dual Taping Station

Model 26

Height choices: 36" or 42".

1" thick high density foam with your choice of optional Bailey upholstery colors, and either Natural or Walnut wood finishes.

Full depth adjustable interior shelf and triple drawers. Spacious 28 quart waste bin, standard Options include Locks and special laminated tops.

Taping Top Dimensions: 26" wide x 32" deep.

Cabinet footprint: 30" deep x 66" wide.

Finished all around; use freestanding as an island.

Weight capacity 500 lbs. per seat.



Model 26

Two units can form a 2 + 2 island!

Model 27 may be combined with or without Model 26 to make any combination of space saving taping stations!



Model 27

Single Taping Station

Model 27

Height choices: 36" or 42",

1" thick high density foam with your choice of optional Bailey upholstery colors, and either Natural or Walnut wood finishes. 32" Deep x 26" Wide Taping Tops; Standard 28 quart Waste bin Supply Top.

Standard Dual adjustable interior shelves & Triple Drawers.

Optional Locks for Doors & Drawers.

Base size 66" Wide x 30" Deep (Model 26 Dual) and 42" Wide x 30" Deep (Model 27 Single).

Finished All-Around, use as an Island or Incline.

Weight capacity 500 lbs. per seat.

Ships Fully Assembled.

Shoulder Wheel Model 605

Ball bearings and brake shoe linings make the resistance mechanism more sensitive to motion. Shoulder wheel is mounted on hardwood frame. Arc of motion ranges from 10" to 39" by moving the handle.

Adjustable in height as much as 26".

All metal parts are chrome plated.

Diameter 37 1/2".

Wall mounted



Model 605



Model 725

Interlocking Stack Stool Set Model 725

Interlocking Stack Stool Features:
Set of Three (3) Stools, 12" W x 24" L
Non-Slip Tread with Natural Wood Construction & Finish.

Three Heights: 4", 6" and 8". Stack Height: 18"
Easy carry handle recess.

Ideal for Step-Ups or as Mini Stairs

Weight capacity 500 lb.

Ship weight 30 lb.

Custom sizes available

Interlocking Stack Stool Set Model 725P

Set of Three (3) Stools, Plyometric 24" x 24" Square, unique design locks together when stacked to prevent lateral movement in use

Non-Slip Tread with Natural Wood Construction & Finish.

Three Heights: 8", 10" and 12", Total Stack Height: 30"

Easy carry handle recess.

Ideal for Step-Ups or as Mini Stairs

Weight capacity 500 lb.

Ship weight 80 lb.



Model 725P

Quadriceps Board Model 761 / 761P

Foam padded and covered with vinyl. (Model 761P)
Standard color is Camel.
Quality hardwood natural finish.

5 1/4" wide x 27 1/2" long x 2" thick when folded.

Model 761 Model 761P



Quadriceps Board Model 7630P / 7640P

Model 7630P

10" wide x 28 3/4" long.
1" padded top with your choice of optional Bailey upholstery colors

Model 7640P

20" wide x 28 3/4" long.
1" padded top with your choice of optional Bailey upholstery colors



Model 7630P

Model 7640P

Adjustable Incline Board Model 7520

14" x 18" foot board with non-slip tread.
Positive lock adjustment handles.
5° - 10° - 15° - 20° - 25° adjustments.
Non-slip bottom surface.
Natural finish.



Model 7520

Ankle Exercise Model 212

This ankle exerciser allows two range of motion exercises. Also promotes safe and proper muscle stretches
12" square x 3-3/4" high
Non-slip safety tread



Model 212

Leg Elevation Board Model 7500

15" wide x 36" long.
7 positions to a 30° incline.
3 pillows with Velcro.
Easy carry hand slot.
Natural finish
Your choice of optional Bailey upholstery colors



Model 7500

Progressive Incline Board Set Model 7540

16 1/4" wide x 60" long.
4 units, each 16 1/4" wide x 15" long with non-slip treads.
1 each 15° - 20° - 25° - 30° incline.
Natural finish


Individual Progressive Incline Boards

Model 7541 (15°)

Model 7541 (20°)

Model 7541 (25°)

Model 7541 (30°)



Model 7540

Crutch and Cane Rack Model 790

Double rods can store up to 30 pair of crutches.
Single rods below are designed to store up to 10 canes. This space saving unit makes it easy to be neat and organized.
Made from ash hardwood
Size 60" long x 17" high and extends, with hangers, 15 1/2" from wall
Mounting hardware not supplied



Model 790

**All Purpose Stools
Models 747, 748 and 749**

Model 747

Height adjusts from 22" to 29"
Requires only 22" of floor space
Seat is 2" thick x 14" diameter
2" ball bearing casters
Standard color is black



Model 747

Model 748 - Stool with 5 Legs

Height adjusts from 22" to 29"
Requires 24" of floor space
Seat is 2" thick x 14" in diameter
5 point base
2" ball bearing casters
Standard color is black



Model 748

Model 749 - With 5 Legs (Not Shown)

Height adjusts from 22" to 29"
Requires 24" of floor space
Seat is 2" thick x 13" square
Back is 10" x 7" x 1-1/2"
Back tilts from 0° to 30°
5 point base
2" ball bearing casters



Model 749

**Step Stool
Model 720**

14" wide x 24" long x 6" high
Polished trim
Stores easily under most treatment tables
Natural finish is standard
Walnut finish is available



Model 720

**Step Stool
Model 721**

12" square x 9" high
Chrome plated, tubular steel frame
Polished trim
Non-skid rubber matting

**P.T. Stools
Model 750 / 752**

Model 750

Chrome plated tubular steel frame
Height 17"

Model 752

Casters are replaced by leg caps
Height 17"



Model 752



Model 750

**Step Stool
Model 722**

12" wide x 16" long x 9" high
Chrome plated, tubular steel frame
Polished trim
Non-skid rubber matting



Model 722



Model 721

**Bariatric Platform Stool
Model 4525**

22" x 22" x 6" high.
Non-slip rubber tread surface.
Dual hand cutouts for portability.
Special Size stores underneath Bariatric Treatment Tables.
Special design allows 1000 lb patient capacity!
Ships fully assembled.
Ship weight 20 lbs.



Model 4525



Whirlpool Chairs

Mid-Size Whirlpool Chair Without Casters Model 7701

Base designed to fit tank contour
Safety Belt
Height adjusts from 23" to 28-1/2"
Seat 18" x 20", back 10" x 16"



Model 7701



Model 7702



Model 7712

Low-Boy Whirlpool Chair Model 7702

Four 3" swivel casters, two with locks
Height adjusts from 22-1/2" to 28-1/2"
Seat 18" x 20", back 10" x 16"

Low-Boy Whirlpool Chair Without Casters Model 7712

Height adjusts from 19-1/2" to 25-1/2"
Seat 18" x 20", back 10" x 16"

Whirlpool Chair Model 7700

Four 3" swivel casters, two with locks
Safety belt
Height adjusts from 33" to 46-1/2"
14" x 19" non-slip rubber foot platform
Seat 18" x 20", back 10" x 16"



Model 7700



Model 7710

Whirlpool Chair Model 7710

Safety belt
Height adjusts from 29-1/2" to 43-1/2"
14" x 19" non-slip rubber foot platform
Seat 18" x 20", back 10" x 16"



Need Whirlpool Tables?
see page 5 for more details

Benches

These multi-purpose benches are a great addition to weight, locker, and training rooms, as well as P.T. and O.T. departments.

Solid ash and steel construction.
Protective floor pads.
Black enamel steel frame.
Natural finish.

Model 3030 - (front bench)

10" wide x 30" long x 18" high

Model 3031 - (middle bench)

10" wide x 48" long x 18" high

Model 3032 - (rear bench)

10" wide x 60" long x 18" high



Stall Bars Model 775

The top rungs are designed for chin-ups and pull-ups. Stall Bars are constructed of the finest hardwood with oval shaped rungs. Sections may be combined to form as large a unit as required.


Width 36", height 96".

Natural wood finish.

Model 3020 - Stall Bar Bench

Bench size 14" x 23" long x 19" high.

Black vinyl is standard with your choice of optional Bailey upholstery colors
Natural finish is standard.
Walnut finish is available.



Model 775

Model 3020

Cuff Weight Carts

Model 783 - Small Cuff Weight Cart

27" wide x 25-1/2" deep x 44" high

Weight capacity 200 lbs.

Natural finish is standard
(cuff weights not included)

Model 784 - Large Cuff Weight Cart

42" wide x 25-1/2" deep x 54-3/8" high

Weight capacity 400 lbs

Natural finish is standard
(cuff weights not included)

Professional Cuff Weight Set Model 1895

Set includes 1 each of the following:

Model 783 cart with 1/2, 1, 1-1/2, 2, 2-1/2, 3, 4, 5, 7-1/2, 10, 12-1/2 and 15 lb. Cuff Weights



Model 783

Model 784

Model 1895





Model 780
Weights Not Included

Model 1890
Weights Included

Professional Fitness Weight Wagon Set Model 1890

This complete set combines the mobility and strength of our Model 780 Weight Wagon with the comfort and color of all ten pair of our Fitbells

Weight Wagon Model 780

All wood mobile storage wagon with chrome plated steel handles and four heavy duty 4" swivel casters. Special slotted top rack stores up to twenty dumbbells, while the four section lower shelf will hold cuff and disc weights.

25" wide x 42" long x 46" high

Solid hardwood construction with multi-layered birch dividers and shelving. Weights Not Included

Cuff Weights

Our Cuff Weights come with a long Velcro® strap to provide a comfortable and secure fit for the wrist or ankle. The Cuff Weights are also color coded since there are 17 weights available and each one has the weight printed, in both pounds and kilograms, on the front.



- Model 1866 - 1/2lb. Walnut
- Model 1868 - 1 lb. Blue
- Model 1869 - 1-1/2 lb. Olive
- Model 1870 - 2 lb. White
- Model 1871 - 2-1/2 lb. Red
- Model 1872 - 3 lb. Gold

- Model 1873 - 4lb. Turquoise
- Model 1874 - 5 lb. Black
- Model 1875 - 6 lb. Beige
- Model 1876 - 7 lb. Lemon
- Model 1877 - 7-1/2 lb. Orange
- Model 1878 - 8 lb. Red

- Model 1879 - 9 lb. Parchment
- Model 1880 - 10 lb. Brown
- Model 1881 - 12-1/2 lb. Olive
- Model 1882 - 15 lb. Tan
- Model 1883 - 20 lb. Blue

Adult Posture Mirror with Floor Stand and Casters Model 700

Glass size 23" x 67"
Overall size 28" x 75" including casters

Model 700



Model 787 Combo Mirror Weight Cart

Footprint 25" x 31" x 80" high, with 4" HD Swivel Casters
Holds standard size products: Six 50 Yard Rolls Exercise Band, Ten Pairs Fitbells, and Twelve Pairs Cuff Weights
Full Size ANSI Z97.1 Safety Glass 28" x 72", Adult Posture Mirror
Balanced and robust construction to prevent tipping
Weight capacity 400 lb.
Ships assembled except for weight hooks; Ship weight 280 lb.
Weights and Bands not included, available separately



Model 787

Fitbells

These solid cast iron dumbbells are provided with a vinyl coating for maximum comfort, safety, and style. (Sold individually)

- Model 1801 - 1 lb. Pink
- Model 1802 - 2 lb. Lime
- Model 1803 - 3 lb. Grape
- Model 1804 - 4 lb. Yellow
- Model 1805 - 5 lb. Blue
- Model 1806 - 6 lb. Orange
- Model 1807 - 7 lb. Dark Green
- Model 1808 - 8 lb. Red
- Model 1809 - 9 lb. Black
- Model 1810 - 10 lb. Royal Blue



Half Roll Positioning Pillows

Model 80 - Half Roll
8" wide x 20" long x 4" high

Model 81 - Half Roll
6" wide x 20" long x 3" high

Model 82 - Half Roll
12" wide x 24" long x 6" high

Model 83 - Half Roll
16" wide x 24" long x 8" high

Model 84 - Half Roll
24" wide x 36" long x 12" high



Multi Set Positioning Pillows
Model 76

One bolster 3" diameter x 12" long
One bolster 6" diameter x 12" long
One rectangle 8" wide x 12" long x 3" high




Bolsters

Model 73 - Bolster
6" diameter x 24" long

Model 74 - Bolster
8" diameter x 18" long

Model 75 - Bolster
8" diameter x 24" long

Model 86 - Triangle Positioning Bolster
15" wide x 16" long x 10" high



Wedge Positioning Pillows

Model 950 - Wedge
7" wide x 17" long x 5" high


Model 951 - Wedge
8" wide x 18" long x 8" high

Model 952 - Wedge
14" wide x 18" long x 4" high

Model 953 - Wedge
26" wide x 24" long x 7-1/2" high

Model 954 - Wedge
26" wide x 24" long x 10-1/2" high

Model 955 - Wedge
26" wide x 24" long x 12-1/2" high



CAN DO®

We carry a full line of Can-Do Exercise products, ask us for details.

Standard Fabric Color Selector - Soft Vinyl



Note: These colors are provided as a guideline and will vary for every printer and monitor. If Color matching is required please request a sample swatch before ordering!

NOTE:

These colors are provided as a guideline and will vary for every printer and monitor.

If color matching is required please request a sample swatch before ordering!

Sherwin-Williams® - Custom Color Wood Matching

Note: We can custom match all your existing or planned stain color by providing us with an actual finished wood sample. Additional costs are associated with this service.

Stain Colors

Wood Classics® Interior Oil Stain is available in 42 colors. Tintable in a wide array of colors for ultimate design flexibility and color matching

Please visit your neighborhood Sherwin-Williams Store.

Limited Warranty

Bailey Manufacturing Company warrants the products we manufacture to be free of defects in materials, workmanship or design under normal use and service conditions for five years after the date of original purchase (two years for Bailey Basics). Resale products and components manufactured by others will carry the warranty of the manufacturer. If the product and/or equipment should become defective within the designated warranty period, Bailey Manufacturing Company will repair or replace it (at our option) free of charge, including return transportation to you, provided you deliver the equipment or product prepaid to Bailey Manufacturing Company at 118 Lee Street, Lodi, Ohio 44254. Please receive authorization before returning product to the factory as Bailey Manufacturing Company will not accept collect shipments. Bailey Manufacturing Company reserves the right to repair or replace component parts without return of the entire unit.



This warranty does not include damage resulting from accident, abuse or misuse of the product and expressly excludes normal wearing of parts or defect caused by transportation, accident, fire, flood, alteration, or negligence.

Bailey's warranty liability is limited to that stated above; and other than the warranty stated above, there are no warranties expressed or implied. Bailey will have no responsibility for consequential or incidental damages.



BAILEY

Bailey Manufacturing Company
P.O. Box 130
Lodi, Ohio 44254-0130



Phone: 1-800-321-8372

Fax: 1-800-224-5390

www.baileymfg.com

BAILEY

Manufacturing Company

The Quality Leader in Sports Medicine Products

