

BAILEY

Model 494 Space Saver Cabinet Treatment Table

Model 494 Features:

- 30" height, 24" wide and 72" long (extended), 45" long folded
- 1" thick vinyl upholstery in US 387 Camel standard, with your choice of optional Bailey upholstery colors, and either Natural or Dark Walnut wood finishes Standard, or any Sherwin Williams color Optional..
- One Drawer and Pass Through Armboard, Powder Coated Stirrups.
- Sliding Door Storage Cabinet
- Weight capacity 300 lb. (center section only)
- Ship weight 200 lb.
- Ships fully Assembled



***Model 494
Shown with Camel Upholstery and
Natural Wood Finish***