

BAILEY

Model 815 ADA Convertible Training Stairs

ADA Staircase Features:

- 36" wide access, with 36" x 36" Platform at 18"H.
- 500 LB capacity, with solid hardwood and engineered hardwood construction.
- Dual handrail design at 35" and 18" for adults and children.
- Step height 6" (3 each) and 4.5" (4 each), tread depth 11", with anti-slip tread and platform surfaces
- Overall Footprint 75" x 86" (angle); 125" x 39" (straight)
- Max height at upper platform handrail is 53".
- Minimum required ceiling height 95" to 101".
- Ships "knocked down", assembly required
- Ship weight 650 lbs.



ADA Training Stairs
*36" Wide - Convertible
Straight Configuration*



ADA Training Stairs
*36" Wide - Convertible
Angle Configuration*