

BAILEY


Model 487 Classroom Treatment Table w/ Removable Mat



Shown w/ Mat



Shown w/o Mat



University & College Health Sciences **Classroom Health Training Equipment**

- Model 487 Classroom Version H-Brace Treatment Table.
- Solid Hardwood and Birch Plywood construction, Natural Finish.
- "English Oak" Laminated Top w/ radius corners.
- Removable Mat w/ Hook & Loop attachment, 2" thick vinyl covered mat with high-density foam, Imperial Blue standard color.
- Table size 30"W x 78"L x 30"H; optional sizes, vinyl colors, and Sherwin Williams wood stain colors are available, ask for current custom pricing!
- Requires on site assembly.

BAILEY

Model 487 Classroom Treatment Table w/ Removable Mat

*Table shown below with mat installed.
Standard floor to table clearance at chair is 22".
Other heights/dimensions available at extra cost.*



Model 487 Classroom Treatment Table w/ Removable Mat

*Table shown below WITHOUT mat installed.
English Oak Laminate is standard.
Other laminates available at extra cost.
Weight capacity 350 LBS.*

