

# BAILEY


## Model 34 / 34A Stirrups

### **Model 34 & 34A Features:**

- Factory Installed, ships assembled or ready to install on tables
- Two Stirrups with 7 positions each; angle does not adjust
- Black Powder Coated finish is non-allergenic, and easy to disinfect and clean
- Fits most Models 400 series, 410 series, 420 series, 430 series, 440 series, and 490 series Space Saver tables.
- NOT for universal or generic fit; may only be used pre-installed on a Bailey table, or as a replacement part for existing Bailey tables.



**Model 34A**  
Replacement Table Stirrups  
w/ Mounting Blocks



**Model 34**  
Optional  
Table Stirrups

**Model 34**  
(installed on table)