


# BAILEY

## Model 355H Horseshoe Activity Table

*Classroom - Rehab - Long Term Care*

### **Features:**

- Multiple places (approx. 4), adjustable height, and wheelchair accessible.
- Height adjusts 22" to 34" with Locking Knobs.
- Overall width 72", overall depth 48", cutout is 36" wide x 30" deep.
- Table surface is uniform 18" wide all around cutout.
- Laminated "English Oak" top with edge banding and rounded corners.
- Reinforced Top at all 5 legs and stress points.
- Weight capacity 150 lb.
- Some assembly required.



See our full catalog of products at  
[www.baileymfg.com](http://www.baileymfg.com)